

BeSure™ Midstream One Step Home Pregnancy Test

- Highly Sensitive
- Ready when You are
- Results in Seconds
- 99% Accurate

This package contains everything you need to perform 1 pregnancy test at home.

BeSure Midstream One Step Home Pregnancy Test

- **Early detection**, QuickStep Midstream can be used as early as 7 days after conception.
- **Convenient**, test any time of the day.
- **One Easy Step**, place the absorbent end in your urine stream for 10 seconds.
- **Easy to read results**. (See diagram)



*Pregnant




*Not Pregnant

*Refer to enclosed instructions.
Always confirm test result with your Doctor.

For external use only. Store at 2 to 30° C.
Contains monoclonal antibodies to human chorionic gonadotrophin.

Absorbent Panel Result Window Easy Grip Handle



This product has been specially packaged to minimise harm to our environment.

Product of Canada. Distributed by QUANTUM MEDICAL
1040 Dandenong Road, Carnegie Victoria 3163, AUSTRALIA

BATCH: 06173C5 EXPIRY: MAY 2007

Instructions for Use

Please read all information in this leaflet before performing the test.

How to Perform the Test

1. Open the sachet and remove the test.
2. Pull the cap off the applicator to reveal the absorbent panel. (See diagram)
3. Hold absorbent panel in urine stream for 10 seconds.
4. Reapply cap to applicator.
5. Wait for the result to appear.

A Positive result can appear between 30 seconds and 5 minutes.
Always wait the full 5 minutes to be certain of the result.
DO NOT interpret results after 10 minutes.

Reading the Result



Pregnant

TWO lines in the result window indicate a Positive result.



Not Pregnant

ONE line in the result window indicates a Negative result.



Invalid

NO lines in the result window indicate an Invalid test.

Always confirm the result with your Doctor.

Questions & Answers

Please read all information in this leaflet before performing the test.

How does the Test work?

BeSure Midstream Pregnancy Test works by detecting the presence of the pregnancy hormone hCG (human Chorionic Gonadotropin) in your urine. This hormone is present soon after conception and progressively increases in quantity during the first 3 months of pregnancy.

When can I do the Test?

You can test any time of the day. You don't have to use first morning urine. Because BeSure Midstream has such a high sensitivity to hCG you are generally able to detect pregnancy as soon as seven days after conception.

How accurate is the Test?

BeSure Midstream is over 99% accurate in laboratory testing.

Does alcohol or any medications affect the Test?

No, not unless they contain the pregnancy hormone hCG.

What do I do if the Test is Positive?

Always consult your Doctor to confirm the test. This test does not intend to replace your Doctor's diagnosis, it is only designed to show whether or not the pregnancy hormone is present in your body.

What do I do if the Test is Negative?

You are probably not pregnant, but if you don't get your period within a couple of days, you may want to re-test. If the re-test is negative and you still have not got your period, you should consult your Doctor.

What do I do if the Test is Invalid?

The most common cause for an Invalid test is when moisture or urine sample gets into the result window. Using another BeSure Midstream test kit, you should repeat the test as instructed on this leaflet with care to avoid any splashes.

What if the one line is faint and the other bold?

This indicates that a very small amount of hCG has been detected, and you are probably pregnant. You may want to wait a few days and do the test again. Always confirm the result with your Doctor.

How long will the result show in the result window?

The result will show for many months, slowly fading over time.

Who do I contact if I have any other questions?

For further confidential information and advice on using BeSure Midstream, please contact:
QUANTUM MEDICAL, 1040 Dandenong Road, Carnegie
Victoria 3163 AUSTRALIA, Tel.: 03 9572 5111 Fax: 03 9572 5199